

Single Cell RNA Analysis

Biocore Bootcamp



University of
Massachusetts
Medical School

Bioinformatics Core

Overview

1) Data Sets


2) Why single cell?

2) inDrops Technology



3) Dolphinnext pipelines

4) Data Structures

BMDC Data Set




Skin Data Set



Why Single Cell??

Solid Diverse Tissue



Minimal Prep



Isolate Cells




Bulk RNA Sequencing




Single Gene Expression Value

Adapted from T. Nawy, Nature Methods, (2014).


Why Single Cell??



Why Single Cell??




Why Single Cell??



Adapted from T. Nawy, Nature Methods, (2014).

Why Single Cell??


- Not all problems necessitate scRNA-seq
- It is well suited when populations are heterogeneous
- It is a powerful tool for studying intra and inter cell type variations in gene expression
- Useful for unbiased discovery



Adapted from T. Nawy, Nature Methods, (2014).


Single Cell Workflows

Single Cell RNA Sequencing Workflow




Single Cell Workflows



Single Cell RNA Sequencing Workflow




inDrops Technology




inDrops Technology



inDrops Technology



inDrops Technology



inDrops Data Processing

Cell
barcode UMI cDNA (50-bp sequenced)

```
AAATTATGACGATGTGCTTG.....GACTGCAC
CGTTAGATGGCAGGGCCGGG.....CTCATAGT
GACCTACGAGTTAGTTTGTA.....GCTCATAA
GTTAAACGTACCTAGCTGT.....GATTTTCT
ACGTCACCTTTGTGGGGT.....ATAAGCTC
TTGCCGTGGTGTATGGAGG.....CCAGCACC
AGTCCATGTGCGGCAGGTTT.....GTTGGCGT
AAATTATGACGAGTTTGTA.....AGATGGGG
CCAAAGATGTCTCTAGGCT.....GGGACGA
GTTAAACGTACCAAGCTTG.....CAAAGTTC
TTTTGACCACTCGTGAGGG.....TTCCAAGG
ACTGTCCATGCCCTGTGTA.....TGGTACGT
CGTAAACCAATAATCCGGTG.....TTAAACCG
.....
.....
.....
```

(Hundreds of millions of reads)

cDNA alignment to
genome and group
results by cell

```
Cell 1 { TTGCCGTGGTGTGGCGGGA.....CGGTGTTA } DDX51
        { TTGCCGTGGTGTATGGAGG.....CCAGCACC } NOP2
        { TTGCCGTGGTGTCTCAAGT.....AAAATGGC } ACTB
        { ..... }

Cell 2 { CGTTAGATGGCAGGGCCGGG.....CTCATAGT } LBR
        { CGTTAGATGGCACGTTATA.....ACGCGTAC } ODF2
        { CGTTAGATGGCATCGAGATT.....AGCCCTTT } HIF1A
        { ..... }

Cell 3 { AAATTATGACGAGTTTGTA.....GGGAATTA } ACTB
        { AAATTATGACGAGTTTGTA.....AGATGGGG }
        { AAATTATGACGATGTGCTTG.....GACTGCAC } RPS15
        { ..... }

Cell 4 { GTTAAACGTACCTAGCTGT.....GATTTTCT } GTPBP4
        { GTTAAACGTACCGCAGAAGT.....GTTGGCGT } GAPDH
        { GTTAAACGTACCAAGCTTG.....CAAAGTTC } ARL1
        { GTTAAACGTACCTTCCGGTC.....TCCAGTCG }
        { ..... }
        { ..... }
```

(Thousands of cells)

Count unique UMIs
for each gene
in each cell

Create digital
expression matrix

inDrops Data Processing

Cell:		1	2	...	N
<i>GENE</i>	1	1	2		14
<i>GENE</i>	2	4	27		8
<i>GENE</i>	3	0	0		1
	:	:	:		:
	:	:	:		:
	:	:	:		:
<i>GENE</i>	M	6	2		0

Dolphinnext Pipelines

To process the FASTQs that the instrument generates into a digital gene expression matrix involves many steps, which can be run as a continuous pipeline with dolphinnext



Dolphinnext Pipelines

<https://dolphinnext.umassmed.edu/>

Dolphinnext Pipelines : bcl2fastq

1) bcl2fastq

- valid sample barcode




```
@NS500602:550:HM2TFBGX5:1:11101:7989:1090:CAACATTGGCCCAAT:CTGACTTA 1:N:0:TATAGACG  
TTCCCNCTAGACATTTCTGTGCATAGATTTTGGTGTGTTTACATAGTCGTTATTCTG  
AAAAA#EEEEEEEEEEEEEEEEEEEEEEEEEEEE/EAAAAEEEEEEEEEEEEEEEEEA6/
```

Dolphinnext Pipelines : extract valid reads

- 1) bcl2fastq
- 2) extract reads


- valid cell barcode
- valid UMI



```
@NS500602:550:HM2TFBGX5:1:11101:7989:1090:CAAACATTGGCCCAAT:CTGACTTA 1:N:0:TATAGACG
TTCCCNCTAGACATTTCTGTGCATAGATTTTGGTGTGTTTACATAGTCGTTATTCTG
AAAAA#EEEEEEEEEEEEEEEEEEEEEEEEEEEE/EAAEEEEEEEEEEEEEEEEEEA6/
```


Dolphinnext Pipelines : alignment

- 1) bcl2fastq
- 2) extract reads
- 3) align



Dolphinnext Pipelines : ESAT

- 1) bcl2fastq
- 2) extract reads
- 3) align
- 4) ESAT



Data Structures : UMI Table

	Cell 1	Cell 2	Cell 3	Cell 4...
Gene 1	42	43	10	9
Gene 2	25	24	2	3
Gene 3	10	9	100	98
Gene 4...	40	39	4	5
SUM	117	115	116	115

Bulk = Average of all genes across all cells

Bulk = proportion of each cell type *
cell types average expression profile

Data Structures : UMI Table

Group A				
	Cell 1	Cell 2	Cell 3	Cell 4...
Gene 1	42	43	10	9
Gene 2	25	24	2	3
Gene 3	10	9	100	98
Gene 4...	40	39	4	5
SUM	117	115	116	115
Group B				

Data Structures : UMI Table

+1000 Cells

+10000 Genes

	Cell 1	Cell 2	Cell 3	Cell 4...
Gene 1	42	43	10	9
Gene 2	25	24	2	3
Gene 3	10	9	100	98
Gene 4...	40	39	4	5
SUM	117	115	116	115

**How can we do this
process in a high
throughput
manner?**

Data Structures : UMI Table

+1000 Cells

+10000 Genes

	Cell 1	Cell 2	Cell 3	Cell 4...
Gene 1	42	43	10	9
Gene 2	25	24	2	3
Gene 3	10	9	100	98
Gene 4...	40	39	4	5
SUM	117	115	116	115

What if we want to store information about these cells or genes??

Data Structures : Expression Set Class

ROW = CELLS
COLS = CELL
METADATA

phenoData

ROW = GENES
COLS = GENE
METADATA

featureData

	Cell: 1	2	...	N
<i>GENE 1</i>	1	2		14
<i>GENE 2</i>	4	27		8
<i>GENE 3</i>	0	0		1
\vdots	\vdots	\vdots		\vdots
<i>GENE M</i>	6	2		0

assay(s)

e.g. 'exprs'

DGE

UMI TABLE

Data Structures : exprs()

Genes	Exprs		
	0hrA_TGACGGACAAGTAATC	0hrA_ATGGGCACACCTTGCC	0hrA_TCGAAGCTGTTGCACG
	0610007P14Rik	0	2.10985244161277
	0610009B22Rik	0	1.72247211812165
	0610009O20Rik	0	0
	0610010B08Rik	0	0
	0610010F05Rik	0	0
Cells			

Data Structures : pData()

pData

	size_factor	UMI_sum	x	y	iPC_Comp1	iPC_Comp2...	iPC_Comp12	Cluster	Timepoint
0hrA_TGACGGACAAGTAATC	1.47	1604.45	0.31	0.72	-0.02	0.01	-0.01	Cluster4	0hr
0hrA_ATGGGCACACCTTGCC	1.21	1333.71	0.93	0.33	-0.02	0.01	-0.01	Cluster3	0hr
0hrA_TCGAAGCTGTTGCACG	1.30	1102.09	0.68	0.13	-0.02	0.02	-0.00	Cluster2	0hr
0hrA_TGTTTGAGTCGGTTCG	1.63	1210.53	0.86	0.32	-0.02	0.02	-0.00	Cluster3	0hr
0hrA_TAAATAGGCACAAGGC	0.43	1714.06	0.90	0.29	-0.02	0.01	-0.01	Cluster3	0hr
0hrA_GATTAGACGGGAACCT	0.64	946.06	0.48	0.59	-0.02	0.01	0.00	Cluster1	0hr

Normalization

tSNE

PCA

Clustering Metadata

Data Structures : fData()

fData

	C1_score	C2_score	C3_score	C1_bulk	C2_bulk	C3_bulk	0hr_score	1hr_score	4hr_score
0610007P14Rik	9704	704	2572	0.05	0.23	0.19	1187	5052	10278
0610009B22Rik	5293	642	1181	0.04	0.11	0.10	503	11045	6766
0610009O20Rik	7535	2732	7733	0.00	0.01	0.01	6310	7579	4016
0610010B08Rik	3184	5176	9845	0.00	0.01	0.00	8161	5170	3772
0610010F05Rik	75	11158	8698	0.11	0.01	0.04	10373	662	4595
0610010K14Rik	3888	5066	2634	0.07	0.10	0.12	1774	4074	10087

Cluster Markers

Aggregated Bulk

Timepoint Markers